



PHILIPS
RESPIRONICS

Sleep therapy

PAP accessories

Rediscover dreams




PAP accessories

for your comfort and freedom to move

Discover new ways to enhance your therapy experience:

- Enjoy a better quality of life with a variety of accessories that allow you to conveniently travel with your sleep therapy system or store components discreetly at home





PAP battery kit

Power your sleep therapy device at home or on the go with the Philips Respironics PAP Lithium ion battery kit. It's designed to provide an uninterruptible power supply (UPS) to your DreamStation or System One sleep therapy devices.

Designed for peace of mind

Using the battery pack's integrated uninterruptible power supply (UPS) will help reduce therapy disruptions when the power goes out. Continuously connecting the battery kit to your Philips Respironics sleep therapy device's power supply keeps it charged and ready for travel or UPS protection.

Convenient for travel

The battery pack is our lightest and smallest ever, offering easy portability so you can continue therapy away from home. It's great for camping, hotel use, etc. The battery pack uses the power supply from your Philips Respironics sleep therapy device to reduce the number of cords needed for travel. And it meets FAA requirements for airplane travel.

Main benefits

- The battery pack provides more than 14 hours of battery life* for PAP device.
- Recharge time: 3-4 hours.
- Charge/discharge indication via the push-button illumination
- It's compact and lightweight at just 1.02 kg.
- It meets FAA requirements for airplane travel (90W lithium ion battery pack).
- No prescription is required.
- It comes with a six-month warranty.
- Its lithium ion battery is environmentally conscious.



Carrying case

Device cable and adaptors

Components

The battery kit accommodates the Philips Respironics DreamStation and System One sleep therapy devices. It includes a rechargeable lithium ion battery pack, a PAP device cable, 2 adaptors (for use with System One 50 series only), and a battery case.

Supports DreamStation, System One 60 series, System One 50 series CPAPs and BiPAPs. Utilizes DreamStation or System One power supply. Additional power supply not included.

Specifications

Dimensions: 1.25" (3.175 cm) x 5.25" (13.335 cm) x 9.0" (22.86 cm)
Weight: 2.25 lbs. / 1.02 kg

Ordering information

PAP Lithium Ion Battery Kit**

Part number

1126942



*PAP run-time test conditions: CPAP mode, pressure 10 cm H₂O, no heated humidification or heated tubing.

**PAP Lithium Ion Battery Kit for use with the Philips Respironics DreamStation and System One CPAPs, BiPAPs, autoSV, BiPAP AVAPS, and BiPAP S/T sleep therapy devices. Includes a rechargeable lithium ion battery pack, a PAP device cable, and a battery case.



DreamStation travel case

The DreamStation travel case* gives you a convenient, compact, and discreet way to travel with your DreamStation sleep therapy device and supplies.

Convenient

The case includes a DreamStation 65W power supply that is lighter and smaller than the standard DreamStation 80W power supply. The 65W power supply powers the DreamStation** device, and the included USB port can be used to charge items such as a smartphone and other electronic devices during travel.

Compact

The lightweight, hard-fabric case offers an ideal combination of protection and portability. Its individual compartments accommodate the DreamStation** device and neatly organize the supplies. An included mesh pouch holds the tubing securely in place.

Discreet

The case is designed to make traveling with sleep therapy equipment less cumbersome. A molded insert liner helps you to thoroughly pack the PAP device and accessories. The case can be attached to a roller suitcase and its zip top facilitates airport security screening. On a plane, it easily fits in the overhead bin or under the seat.

Main benefits

- The hard-fabric case provides added peace of mind that the sleep therapy device will be protected during travel.
- The included power supply offers travel convenience.
- The USB port on the 65W power supply reduces the need for an additional charger.
- The molded insert liner is designed to help you thoroughly pack the PAP device components.



Components

The case accommodates your Philips Respironics DreamStation sleep therapy device (blower only, no humidifier), your mask and tubing, a 65W power supply w/ USB port, and a removable luggage strap.

Specifications

Dimensions: 10.42" (26.46 cm) x 3.94" (10 cm) x 13.78" (35 cm)

Ordering information

DreamStation Travel Kit
(incl. 65W power supply w/ USB port)

Part number

EU1120135

*DreamStation, mask, and tubing not included
**Without humidifier



PAP travel briefcase with detachable laptop bag

The PAP travel briefcase* gives you an all-in-one carrying option that accommodates most PAP therapy systems and humidifiers, a laptop, and other small items. It makes it easy and efficient to travel with a sleep therapy system.


The travel briefcase features two bags that conveniently zip together into one carry-on. The therapy device and accessories fit snugly into customizable compartments on one side of the case, while the laptop and its accessories fit securely into the ample-sized pockets of the other side. There's even space for books, magazines, and other small travel items.

Main benefits

- The case's compartments protect the sleep therapy device and accessories during travel.
- Its zip top facilitates airport security screening.
- The 2 bags zip and unzip easily so either bag can be carried separately.
- It fits within the FAA maximum carry-on size.
- The lightweight case offers easy portability.
- It's designed to attach to a roller suitcase.




Power supply



Tubing



DreamStation with humidifier**



Removable laptop bag



Detachable shoulder strap

Components

The briefcase accommodates most PAP therapy devices and accessories and standard-size laptop computers.

Specifications

Dimensions: 16.14" (40.99 cm) x 11.5" (30.5 cm) x 7.5" (19.99 cm)
Color: Black

Ordering information

PAP/Laptop travel briefcase

Part number

1114784

*DreamStation, power supply, and tubing not included.

**Or other PAP therapy system with humidifier.



Mask and tubing bedside organizer

The bedside organizer gives you a way to conveniently store your mask and tubing without cluttering your nightstand – both at home and during travel. It also keeps the components dust-free and ready for use.

The organizer has a weighted flap to allow it to hang over a nightstand or in between mattress and box spring. It even features a zippered compartment for the mask and tubing, making the organizer a great travel companion.

Main benefits

- The organizer provides a clean, discreet place to store the mask and tubing while not in use.
- Its design looks right at home in the bedroom.
- It can be used on a nightstand or in between a mattress and box spring.

Specifications

Dimensions: 7.99" (20.3 cm) x 10.94" (27.8 cm) x 3.74" (9.5 cm)

Color: Light gray

Ordering information

Bedside organizer

Part number

1122135



LiquiCell® nasal CPAP cushions

LiquiCell® is a thin, pliable cushion infused with fluid that can be applied virtually anywhere on the face, making it versatile. The cushion enhances comfort by acting as a barrier between the user's face and the mask. Gentle adhesive keeps the cushion securely in place.

Ordering information

LiquiCell nasal CPAP cushions, pack of 30

Part number

1109097

*LiquiCell is a registered trademark of Birchwood Laboratories, Inc.

